


Since 1951 the name ESCHWEILER has stood for innovative technology in the area of blood gas and electrolyte analysis.

Reliability and simple operation are our goals – the consistent fulfilment of these objectives has distinguished blood analysers from ESCHWEILER for decades.



#### TECHNICAL DATA

measured parameter	range	resolution
pO <sub>2</sub>	0 - 800 mmHg (SI-units selctable)	0,1 mmHg
pCO <sub>2</sub>	5 - 200 mmHg (SI-units selctable)	0,1 mmHg
pH	6,000 - 8,000	0,001 pH
total-hemoglobin (tHb)	4 - 30 g/dl	0,1 g/dl
barometric pressure	100 - 900 mmHg (SI-units selctable)	1,0 mmHg
Na <sup>+</sup>	20 - 250 mmol/l	1,0 mmol/l
K <sup>+</sup>	0 - 20 mmol/l	0,01 mmol/l
Ca <sup>++</sup>	0 - 5,0 mmol/l	0,01 mmol/l
Li <sup>+</sup>	0,4 - 5,0 mmol/l	0,01 mmol/l
Cl <sup>-</sup>	20 - 250 mmol/l	1,0 mmol/l
glucose	0 - 30 mmol/l (0-540 mg/dl)	0,1 mmol/l
lactate	0 - 20 mmol/l (0-180 mg/dl)	0,1 mmol/l
input parameter		
patient temperature	13 - 43 °C	0,1 °C
hemoglobin (tHb)	0 - 30 g/dl (if not measured)	0,1 g/dl
fraction of inspired oxygen (FIO <sub>2</sub> )	15 - 100 %	only relevant for AaDO <sub>2</sub>
respiratory quotient (RQ)	0,7 - 1,0	only relevant for AaDO <sub>2</sub>
calculated parameter		
hydrogen ion conc. (H <sup>+</sup> )	10 - 1000	0,1 nmol/l
actual bicarbonate (HCO <sub>3</sub> -A)	10 - 50	0,1 mmol/l
standard bicarbonate (HCO <sub>3</sub> -S)	10 - 50	0,1 mmol/l
base excess (BE)	-25 - 25	0,1 mmol/l
standard base excess (SBE)	-25 - 25	0,1 mmol/l
total CO <sub>2</sub> (TCO <sub>2</sub> )	10 - 50	0,1 mmol/l
buffer base (BB)	0 - 100	0,1 mmol/l
O <sub>2</sub> saturation of hemoglobin (O <sub>2</sub> sat)	20 - 100	0,1 %
O <sub>2</sub> content or concentration (O <sub>2</sub> CT)	0 - 40	0,1 %
partial O <sub>2</sub> -press. at 50% O <sub>2</sub> -sat (P50)	10 - 50	0,01 mmHg
alveolar to arterial oxygen-tension grade (AaDO <sub>2</sub> )	0 - 800	0,1 mmHg
anion gap (A-GAP)	0 - 99	0,1 mmol/l
SHUNT	0 - 50	0,1 %
acid base status	relevant diagnosis recorded on printer	
hematocrite (Hct)	0 - 100 %	0,1 %
data input		data output
touch screen	display	10.4" color LED-backlight display (TFT 800 x 600)
barcode reader	printer	fast, low-noise thermoprinter
	interfaces	network (LAN), USB, RS 232, HL7
	data base	integrated data bases for patient- and QC-data
calibration		specimen
automatic calibration	90 min	specimen container
economy mode	240 min	specimen material
sleep mode	no calibration	capillary, syringe and others
		whole blood, serum, plasma and respiratory gas
dimensions/weight		electrical data
height	520 mm	mains voltage
width	450 mm	115 V resp. 230 V
depth	415 mm	frequency
weight	appr. 17 kg	50/60 Hz
		ambient temperature
		12 - 32 °C



modular **pro**

BLOODGAS- AND  
ELECTROLYTE ANALYSIS  
RELIABLE  
EASY-TO-USE





Reliable, fast and easy to use

The new model series ESCHWEILER modular pro offers the user the proven ESCHWEILER concept with state of the art technology:

- modular construction with flexible sensor configuration
- up to 11 directly measured parameters:  $pO_2$ ,  $pCO_2$ , pH,  $K^+$ ,  $Na^+$ ,  $Ca^{++}$ ,  $Li^+$ ,  $Cl^-$ , glucose, lactate, tHb
- ESCHWEILER CALPACK calibration cassette system
- larger touch screen with clear operating instructions
- universal sample port for syringe and capillary
- barcode reader for entering sample data
- management of patient and QC data
- network access via LAN (HL7)



#### Touch screen

The large 10.4 inch touch screen enables fast and simple operation of the ESCHWEILER modular pro in daily use.

- simple and clearly-structured menu navigation
- large operation keys for measuring routines
- extensive service and diagnosis programs
- integrated databases for the acquisition and analysis of patient and QC data

#### Sample port

The new universal sample port of the ESCHWEILER modular pro actively draws the sample to be measured.

Whether capillary or syringe, after insertion of the sample container the ESCHWEILER modular pro takes over all further steps of the sample entry.

#### Calibration

ESCHWEILER CALPACK – the new calibration concept for ESCHWEILER modular pro analysers.

All the necessary calibration solutions for blood gas, electrolyte and metabolite sensors are included in an easy to exchange calibration cassette.

- RFID system for automatic transmission of calibration values
- automatic fill level monitoring by the analyser
- simple menu-driven CALPACK exchange
- proven ESCHWEILER liquid calibration of the gas sensors



#### Sensors

For more than six decades the name ESCHWEILER has stood for high quality and reliability in blood gas and electrolyte sensors. Measuring precision and durability of measuring units are the hallmarks of ESCHWEILER sensors. Low maintenance and exchangeable membranes – this keeps service costs low.

11 sensors for direct measurement are available and can be combined individually according to the requirements of the customer. In addition there are 15 further automatically transmitted calculation parameters. (see "Technical Data")

#### Network (LAN)

- Connection of the ESCHWEILER modular pro to internal clinic networks via LAN interface (HL7)
- Transfer of patient and QC data to the central data administration via ODBC (open database connectivity)

#### Barcode Reader

- Read in coded patient data
- Read in coded QC samples

